TEST METHODS MANUAL

Laboratory Services Division Bureau of Soils and Water Management

Department of Agriculture

TM-LSD-04-16

SECTION : SOIL CHEMISTRY	Issue No.: 3	Effective date: July 1, 2022
SUBJECT : SOIL SALINITY/ALKALINITY - SOIL pH (Alkalinity)	Revision No.: 0	Page 1 of 4

SCOPE

This method is for the measurement of pH on soils that contain excessive concentrations of either soluble salts or exchangeable sodium or both.

PRINCIPLE

Soil pH value is determined with an H^+ ion-selective glass electrode (referenced with a standard calomel electrode) inserted into the soil filtrate from the saturated soil extract. The difference between the H^+ ion activities in the soil solution and the glass electrode gives rise to an electrometric potential difference that is related to the soil solution pH, or the active acidity.

It is determined by the concentration of hydrogen ions $[H^+]$. It is measured using the formula of negative logarithmic scale of the hydrogen ion concentration $[H^+]$, pH = $-\log[H^+]$. The soil pH value is whether an acidic solution (less than 7) or an alkaline solution (greater than 7). So, as the hydrogen ion concentration goes up, the pH value goes down, which makes it more acidic, and vice versa with alkaline solution.

Soil pH is influenced by both acid and base-forming cations (positively charged dissolved ions) in the soil. Common acid-forming cations are hydrogen (H^+), aluminum (Al^{3+}), and iron (Fe^{2+} or Fe^{3+}), whereas common base-forming cations include calcium (Ca^{2+}), magnesium (Mg^{2+}), potassium (K^+) and sodium (Na^+).

The pH of soil changes over time, by monitoring and comparing the past and present soil test there is a probability to see if the soil acidity is increasing over time. Thus, we can prevent this trend from continuing.

TEST PRECAUTIONS

The pH values are dependent on temperature. It should be measured at constant temperature. If a precision greater than 0.1 pH unit is desired, the temperature of the standard solutions, the glass electrode, and the test solution must be within 2° C of one another.

Keep electrodes wet by returning them to storage solution whenever pH meter is not in use. Recommended solutions for short-term storage of electrodes vary with the type of electrode and manufacturer.

Clean the glass electrode periodically if necessary, by soaking it in a solution of a mild detergent and wiping it with a piece of tissue paper, to remove the film of soil that adheres to the electrode surface. A film of $CaCO_3$ that sometimes forms on the surface of the glass electrode during measurement of pH of calcareous soils may be removed by dipping the electrode for a few seconds in dilute HCl and rinsing it with water.

Replace electrodes when pH calibration buffers exceed ± 0.05 pH units, even after recalibration, or when pH of quality control samples consistently exceeded expected error.

Prepared by:	Reviewed by:	Approved by:		
embgamboa	JP Sanche	AW .		
EZRA MÁÉ B. GAMBOA	FLORFINA P. SANCHEZ	GINĂ PUNILO, Ph.D.		
Document Controller	Head, Soil Chemistry Section	Quality Manager		
THIS IS A CONTROLLED DOCUMENT MAINTAINED ELECTRONICALLY.				
THIS IS UNCONTROLLED WHEN DOWNLOADED, PRINTED AND PHOTOCOPIED.				

- **Limited Preview**

For full access, kindly fill out the Document Request

- Form.

Google Form link:

HEALTH AND SAFE https://forms.gle/RbCgCdA54prTS6oN7

Thank you!