### TEST METHODS MANUAL

Laboratory Services Division Bureau of Soils and Water Management Department of Agriculture

## TM-LSD-04-20

SECTION : SOIL CHEMISTRY	Issue No.: 3	Effective date: July 1, 2022
SUBJECT : SOIL SALINITY/ALKALINITY - CHLORIDES (SILVER	Revision No.: 0	Page 1 of 4
NITRATE METHOD)		

#### SCOPE

This method is for the determination of chloride on soils that contain excessive concentrations of either soluble salts or exchangeable sodium or both

#### PRINCIPLE

Chloride in water and in soil can be determined by Silver Nitrate Method. It is based on the principle that in neutral or slightly alkaline sample, potassium chromate can indicate the endpoint of the silver nitrate titration of chloride. Silver Chloride is quantitatively precipitated before red silver chromate is formed.

#### **TEST PRECAUTIONS**

The end point is sharper if the titration is done under a yellow light.

#### EQUIPMENT

a) Balance, precision of at least 0.001 g

#### LABORATORY WARE

- a) Erlenmeyer Flasks, 250Ml
- b) Pipettes, graduated
- c) Beakers
- d) Wash bottle
- e) Microburette, 10 mL
- f) Volumetric Flasks, 1000 mL., 100 mL.
- g) Pipettor

#### **CHEMICALS AND REAGENTS**

- a) Potassium chromate, 5% solution. Dissolve five (5) g of  $K_2CrO_4$  in 50 mL distilled water add a solution of AgNO<sub>3</sub> dropwise until a slight permanent red precipitate is produced. Stand overnight and filter. Dilute to 100 mL with deionized water.
- b) Standard Silver Nitrate, 0.005 M. Dissolve 0.8495 g of AgNO<sub>3</sub> in deionized water to exactly one (1) liter. Keep in brown bottle, away from light. To standardize, titrate ten (10) mL of 0.01 M KCl (0.1864 g of KCl dissolved in 250 mL water) against AgNO<sub>3</sub> solution until the appearance of reddish brown precipitate. Use four (4) drops of K<sub>2</sub>CrO<sub>4</sub> as indicator. Compute for concentration (M) of standard AgNO<sub>3</sub> as follows:

Prepared by:	Reviewed by:	Approved by:		
embgamboa	JP Sanchy	(AW		
EZRA MÁÉ B. GAMBOA	FLORFINA P. SANCHEZ	GINĂ PUNILO, Ph.D.		
Document Controller	Head, Soil Chemistry Section	Quality Manager		
THIS IS A CONTROLLED DOCUMENT MAINTAINED ELECTRONICALLY. THIS IS UNCONTROLLED WHEN DOWNLOADED, PRINTED AND PHOTOCOPIED.				

### **TEST METHODS MANUAL**

aboratory Services Division Bureau of Soils and Water Management Department of Agriculture

# TM-LSD-04-20

Molarity (M) of  $AgNO_3 =$ 

M<sub>KCI</sub> x V<sub>KCI</sub> V<sub>AgNO3</sub> used

#### HEALTH AND SAFETY

Wear proper PPE. Use laboratory coat, close shoes, gas masks or dust masks and appropriate gloves when performing chemical analysis to mitigate Limited Preview use on chemicals.

For Laboratory Use: always check the SDS and COA for each delivered chemical as to confirm if it is compliant

## For full access, kindly fill out the Document Request

Observe careful and proper handling of chemical **Form**ing inorganic salts, strong alkali bases, strong acids and oxidizing agents. Potential direct body contact (e.g. ingestion inhalation) may cause severe irritation and inflammation to the skin, eyes, respir **Google Form link:** re exposure to these corrosive substances may cause serious hattps://forms.gle/RbCgCdA54prTS6oN7

Avoid mixing incompatible chemicals to reduce risks of fire and explosion inside the laboratory. Keep locked up and away from incompatibles such as moisture.

 $K_2CrO_4$  may cause an allergic skin reaction, genetic defects, cancer and damage to organs (blood, kidneys, liver) through prolonged or repeated exposure, very toxic to aquatic life with long lasting effects Precautionary statements (GHS-US): P201- Obtain special instructions before use. P202- Do not handle until all safety precautions have been read and understood Thank you! mist. Contaminated work clothing should not be allowed out of the workplace and avoid release to the environment.

#### PROCEDURE

- 1. Prepare a blank.
- 2. Add four drops of  $K_2CrO_4$  indicator, while swirling, to the blank and to the sample preserved from the carbonate-bicarbonate determination.
- 3. Titrate under a bright light with 0.005 M AgNO<sub>3</sub> solution from a 10 mL microburette to the first permanent reddish brown color. The titration blank correction with the volume of the sample at the end point usually increases regularly from about 0.03-0.20 mL as the volume increases from 2-12 mL.

nte: (1) Titrate from yellow until dark orange color (2) Titrate blanks first to establish the point of color of change.